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Instructions for reprocessing of medical devices produced by vitaclinical

vitaclinical is a business division of
VITA Zahnfabrik (following called vitaclinical)

Manufacturer: VITA Zahnfabrik, Spitalgasse 3,
79713 Bad Säckingen, Germany

General requirements

Deviations from these instructions (for manual and mechanical cleaning and disinfection as well as during sterilisation process) are in the responsibility of the users. Vitaclinical does not bear any accountability for defects, caused by any deviation from the following instructions.

vitaclinical reserves the right to make amendments to this direction due to any new findings.

These reprocessing instructions are available as download (German version) within www.vita-zahnfabrik.com or can be requested directly from VITA Zahnfabrik H. Rauter GmbH & Co. KG, Ballyweg 6, 79713 Bad Säckingen, Germany.

Delivery of medical device is in unsterile condition and must be cleaned/disinfected and sterilised before usage. Used devices as well as devices taken out of sterile packaging are rated as unsterile and as consequence are rated as contaminated and must be treated hygienic and reprocessed. Before using a sterile medical device, the sterile packaging has to be checked for damages. If the sterile packaging is damaged, the enclosed device is to rate as not sterile.

Only validated procedures are allowed to use for cleaning/disinfection and sterilisation. In each cycle the validated parameters must be followed and the used devices must be maintained regularly. If other procedures are applied than the procedures described in this direction, these must be validated through the person in charge for reprocessing.

Handling of used, unsterile medical device from point in time of application and until removal of the steriliser has to happen in order that whether patient, neither user nor any third party is at any risk. Staffs have to wear protective clothing during the whole reprocessing procedures.

Please consider as well country specific regulations and laws as well as hygienic rules for hospitals and doctor's practices. This applies particularly to the different rules regarding effective prion deactivation.

Based on the foreseen purpose (medical devices for surgical and paradontological procedures, they may penetrate the skin and mucosa and come into contact with blood or tissue) we recommend to rate the most critical medical device torque.wrench with complex geometry (inner cavity) as risk group "critical B" according RKI-Directions. Concurrently the product torque-wrench is the "worst-case" product out of the whole sets (single component with lumen/cavity).

According recommendation of Robert-Koch-Institute and for better validationability, for cleaning and disinfection always a mechanical process must be used.

For a high standard of certainty (drying of contamination) the pre-treatment is to start immediately, not later than 30 minutes after usage. It is to start with cleaning/disinfection latest 2 hours after usage.

Dismountable products (torque.wrench) are to disassemble before cleaning and disinfection and to be assembled as well as to be checked for operational availability before sterilisation. For the product torque.wrench we recommend a pre-treatment of lumen (hole). Depending on degree of contamination, we recommend in addition an ultrasonic cleaning. All components surgical.tray consists of should be cleaned machinable and disinfected after usage.

1. Pre-treatment

Immediately after usage (latest 30 minutes after usage) pre-treating of the product should be started.

Process:

- Dismount product as far as possible
- All major visible contamination is to clean under cold running water by using nylon brush or fuzz-free material. Pretreatment could be done in a disinfection lotion (if a critical biological contamination exists).
- For devices with inner lumen (torque.wrench) cavity should be flushed 5 times with single use syringe (5 to 10 ml) in addition to mechanical treatment. Needle should not be made of metal in order to avoid damages at the device.
- Pre-treatment cleaning supplies should be free of aldehyde and alcohol (danger in fixing of proteins).
- Use only protein non-fixing cleaning and pre-treating supplies

Disassembly – for cleaning and disinfecting the torque.wrench:

1. Turn adjuster (auch: torque set screw) with one surface upwards
2. Pull the bolt with spring scale and adjuster out of the handle

Assembly – for sterilizing:

1. Turn adjuster (auch: torque set screw) with one surface downward
2. Insert bolt with spring scale and adjuster into the handle
3. If resistance is felt, turn adjuster on opposing surface and insert

2. Ultrasonic Cleaning (optional)

If strong contamination exists, an ultrasonic cleaning bath could be done in advance. Adhesive contamination on product's surface will be dissolved and could be removed easier, if an ultrasonic cleaning bath is done in advance.

Pay attention to the following when doing ultrasonic cleaning:

- Stated residence time as well concentration given by the producer of the cleaning additive must be followed.
- Cleaning material must suite the device.
- Consider compulsory fluid fill quantity in ultrasonic cleaning bath (e.g. filling up until identification mark).
- Devices must be covered totally by cleaning lotion.
- In order not to affect the ultrasonic effect, devices should only be placed on sieve bowls
- The cleaning lotion must be renewed regularly. Visible contamination is a criterion. In minimum once a day, it must be changed.

Pre-treatment and ultrasonic cleaning will be repeated until no contamination is visible on the device anymore and until the flushing lotion is clear.

3. Mechanical Cleaning/Disinfection (Disinfector/Cleaning- and Disinfection Device RDG) Process:

- Put the devices in disinfector basket respectively in sterilisation sieve (devices should not touch).
- Put disinfector basket respectively sterilisation sieve with devices into disinfector.
- Start the cleaning program.
- At the end of the program, take out disinfector basked respectively sterilisation sieve with the devices. Devices and cavity (eg. torque.wrench) must eventually be dried with filtered compressed air.
- Assemble the devices
- Check and pack the devices immediately, if possible.

Pay attention to the following when choosing a disinfector:

- Disinfector must be maintained and tested regularly.
- Disinfector must follow requirements of DIN EN ISO 15883 and AAMI ST15883 and must have a tested effectivity (e.g. CE-, DGHM/VAH- or FDA-proved)
- Eventually use for thermal disinfection a tested disinfection program (min. A0-Value > 3000 (resp. for old tools until min. 10 minutes at 93°C)). Because with a chemical disinfection the risk exists, that disinfection lotion residuals remain on the device.
- The used program must be suitable for the device and must have a sufficient number of flushing cycles.
- It is allowed to use sterile or low-germ water (max. 10 germs / ml) as well as low-endotoxin water (max. 0,25 endotoxin items / ml, eg. Aqua purificata / Aqua purificata valde (purified water (PW)/highly purified water (HPW))).
- Air used for drying must be filtered respectively must follow requirements given by hospitals.
- Cleaning supplies must be suitable for devices made of plastics and metal.
- While applying a chemical disinfection a suitable disinfection supplies with tested effectiveness (e.g. CE-, DGHM/VAH- or FDA-proved), which is compatible with the used cleaning supplies and the device (refer to material's consistency) should be used.
- Disinfection supplies should be free of aldehyde and alcohol (danger in fixing of blood).
- A sufficient neutralisation according supplier's recommendation is to make if an alkaline cleaning supply was used.
- You must stick to the supplier of the cleaning- and disinfection supplies recommended concentration.

4. Manual Cleaning/Disinfection

(A manual cleaning and disinfection is based on the low effectiveness and worse reproducibility, only to apply if an automated procedure is not available. Pre-treatment should be carried out before manual cleaning/disinfection.)

Manual Cleaning:

If no pre-cleaning was done, device must be disassembled before starting manual cleaning. Hard to clean devices (e.g. inner lumen) should be pre-cleaned before starting manual cleaning.

Process:

- Put the devices in cleaning lotion. Follow residence time and concentration instructions given by the producer of the cleaning supplies.
- Devices with lumen must be flushed at the beginning and at the end of residence time in min. five times with 5 to 10 ml cleaning lotion. To do so, use a clean, unused one-way gun, if possible with a plastic needle.
- Visible, adhesive contamination will be removed while using a soft nylon brush or a soft, fuzz-free cloth.
- Flush devices in min. five times intensive in freshly distilled or deionised water (flush each lumen in addition five times with 5 to 10 ml).

Cleaning process will be repeated until no contamination is visible on the device anymore and until the flushing lotion is clear.

Consider the following points when choosing cleaning supplies:

- Cleaning supplies must be suitable for cleaning of the devices and must be compatible with the devices (suitable for cleaning of metal and plastics). Due to the eventually occurring strong contamination, usage of combined cleaning-/disinfection supplies is not recommended.
- Cleaning supplies should be suitable as well for the ultrasonic cleaning (no foaming while applying ultrasonic).
- Cleaning supplies must have a tested efficiency (e.g. CE-, DGHM/VAH- or FDA-proved).
- You must stick to the supplier of the cleaning- and disinfection supplies recommended concentration and residence time. Only lotions produced freshly are allowed to be used.
- For all flushing steps, only use sterile or low-germ water (max 10 germ / ml); consider a very low endotoxin- (max 0,25 endotoxin items / ml) and particle pollution (e.g. Aqua purificata/ Aqua purificata valde (purified water (PW)/highly purified water (HPW)) according to demand of Pharm. Eur. or USP)).

5. Manual Disinfection

Process:

- Put the devices in disinfection lotion. Follow residence time given by the producer of the disinfection lotion.
- Devices with lumen must be flushed at the beginning and at the end of residence time in min. five times with 5 to 10 ml disinfection lotion. To do so, use a clean, unused one-way gun, if possible with a plastic needle.
- Flush devices in min. five times intensive in freshly distilled or deionised water (flush each lumen in addition five times with 5 to 10 ml).
- Dry devices in a clean, dust-free place by using filtered compressed air.
- Assemble the devices.
- Check the functionality of devices as well as damages occurred. (refer to article „Testing“)
- Pack devices and prepare for sterilisation.

Consider the following points when choosing disinfection lotion:

- Cleaning supplies must be suitable for disinfection of the devices and must be compatible with the devices (suitable for cleaning of metal and plastics). Due to the eventually occurring strong contamination, usage of combined cleaning-/disinfection supplies is not recommended.
- Disinfection supplies must have a tested efficiency (e.g. CE-, DGHM/VAH- or FDA-proved) and be compatible with cleaning supply.
- You must stick to the supplier of disinfection supplies recommended concentration and residence time. Only lotions produced freshly are allowed to be used.
- For all flushing steps, only use sterile or low-germ water (max 10 germ / ml); consider a very low endotoxin- (max 0,25 endotoxin items / ml) and particle pollution (e.g. Aqua purificata/ Aqua purificata valde (purified water (PW)/highly purified water (HPW)) according to demand of Pharm. Eur. or USP)).

6. Testing

After each clean up / disinfection the devices must be clean and free of contamination when making optical test (free of e.g. visible protein residuals or other contamination. If this is not the case, the whole cleaning / disinfection process, if necessary pre-treating, has to be repeated. Worn out, damaged or corroded devices must be sorted out because functionality is affected or other products could be affected.

Also devices with contamination that cannot be removed have to be sorted out.

7. Re-use

Devices which can be rated as full functional after using, reprocessing and testing, can be reused up to five hundred times.

8. Packaging

All devices separated from each other, cleaned and disinfected, will be sorted into foreshen tray and packed in sterile packaging suitable for steam sterilisation (single or double sterile pouches) or will be packed in sterilisation container.

Packaging for steam sterilisation must be steam permeable and heat resistant at min 137°C (279°F), plus must show after sterilisation an efficient protection against recontamination. They must follow the standard requirements: EN ISO 11607 / ANSI/AAMI ST79 / AAMI TIR12:2010 / DIN EN 868-8.

9. Sterilisation

Only cleaned and disinfected devices are allowed to be sterilised. For sterilisation only the listed sterilisation processes are allowed to be used. Other sterilisation processes must be validated by the user.

The given sterilisation-times and –temperatures are minimum requirements and it is not allowed to underturn them. If sterilisation parameter is underfunded, sterility of the device is not guaranteed any more. Such a sterilisation process must be validated by the user.

Exceeding the indicated sterilization temperatures and times does not influence the sterility of the device negatively. These factors can cause a higher material imposition and as conclusion the device is aging earlier. Application of other sterilisation processes (e.g. hot air-, ethylene oxide-, formaldehyde-, radiation- and low temperature-plasma-sterilisation) is not in the producer's responsibility and should be validated by the user before sterilisation.

After sterilisation is completed, the medical devices are to be stored dry, dust-free and protected from solar irradiation. Consider specifications given by producer of sterile packaging in regards to durability. On sterile packaging a sterilisation indicator (as reference to the sterilisation) as well as a sterilisation date must be existent.

Note:

The tray is no sterile container.

Consider the following points when choosing sterilisation process:

- Fractionate vacuum procedure (with sufficient product dryness). Application of less efficient gravity method must be ensured with an additional validation (if needed longer sterilisation time is necessary). For lumen instruments (e.g. needles, handles) this gravity method is unsuitable.
- Steam sterilization devices equivalent DIN EN 13060 respectively DIN EN 285 and DIN EN (ANSI AAMI) ISO 17665-1 validated (valid consignment and product specific assessment of performance).
- Maximum sterilization temperature 137°C (add. tolerance equal DIN EN (ANSI AAMI) ISO 17665-1.
- Sterilization times:
 - Minimum 20 minutes at 121°C, alternatively
 - Minimum 5 minutes at 132°C / 134°C

(Based on recommendation by Robert Koch-Institute, sterilization time at 134°C is minimum 5 minutes). This recommendation presumes the use of an alkaline cleaner with pH greater than 10. If no alkaline cleaner is used, sterilization time should be extended to 18 minutes.

Validation, Evidence, Tests

Documented evidence of applicability and reproducibility of the processes described above such of efficient mechanical cleaning and disinfection, manual cleaning and disinfection as well as steam sterilization was provided by an independent test laboratory. While testing following devices and instruments were used:

- Disinfectant type G 7883 CD (Miele) TDA, serial number 00/74312631 for cleaning and thermal disinfection, steam sterilization device BBC/Technic Labo serial number 5119.
 - Cleaning and disinfection supplies (manufacturer Schülke & Mayr): GIGASEPT INTRU AF for pre-cleaning and manual disinfection as well GIGAZYME for manual and mechanical cleaning.
- The described processes were considered.

Explanation of symbols used

 REF	Order number		Expiration date
 LOT	Batch designation		Store away from sunlight
 STERILE EO	Sterile packaging, EO sterilised		0124 The CE mark confirms that the requirements of the Council Directive 93/42 EEC concerning Medical Devices are fulfilled
	Do not use in the case of damaged packaging		Attention – please refer to the enclosed documents
	Not for re-use		Manufacturer
	Date of manufacture		