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Principles of treatment planning

The surgical part of implant rehabilitation must be preceded by comprehensive patient evaluation, preoperative diagnostics and treatment planning. Deficits in treatment planning can lead to implant loss.

Possible complications

Transitory pain, swelling, phonetic difficulties or gingival irritations may occur.

The following long-term adverse reactions may occur:

Chronic pain associated with the dental implant, permanent paraesthesia, dysaesthesia, bone resorption of the upper / lower alveolar ridge, localised or systemic infections, oronasal or oronasal fistula, irreparable damage to adjacent teeth, implant, jaw, or bone, dental restoration fractures, aesthetic problems, nerve damage, hyperplasia, implant loss.

Bone preparation prior to insertion of a ceramic implant CI

The recommended size of the drills / thread cutters is specified to help achieve sufficient primary implant stability. Hence the diameter of the various profile.drills is 0.4 mm smaller than the outer diameter of the implant. Because ceramic implant is not designed as a self-tapping implant, the thread cutter must be inserted as the final preparatory instrument.

Insertion of a ceramic implant CI

- In the case of mechanical insertion, set a low rpm of up to 10Upm and a maximum torque of 35Ncm (minimum 25Ncm).
- In the case of manual insertion, first screw into place by hand. As the resistance increases, place the torque.wrench onto the hand.wrench, apply slight pressure to the hand.wrench using your finger and use the torque.wrench to slowly screw into position. Max. torque 35Ncm (minimum 25Ncm).

Important: When the ceramic implant CI is in its final position, the upper thread groove must be sunk into the bone.

Healing phase of ceramic implant CI with cer.face 14®

If used for the recommended indications and given sufficient primary stability, an 8-week healing period is necessary for a lower jaw ceramic implant with cer.face 14. A 16-week healing period is required for upper jaw applications.

A longer healing period is necessary when using bone augmentation procedures or in case of incomplete bone contact with cer.face 14.

Modifications of ceramic implant CI abutment

If required, a slight reduction in height can be carried out on the ceramic implant abutment. We recommend the following procedure:

- A maximum of 1 mm height can be reduced occlusally
- Use only a wet grinding turbine
- Apply only small amounts of pressure
- Fine-grained diamonds (as new) with a red color coding (fine 27 – 76 µm) should be used
- An extra fine (yellow 10 – 36 µm) or ultra fine (white 4 – 14 µm) grain should be used during the final milling step
- After abutment modification an impression must be taken either with set.cap or individually in order to transfer the modified shape reduction into the master cast. The use of insert.transfer in this case leads to an inaccurate master cast.

ceramic implant CI impressions

In cases where an impression is planned directly following surgery, insert.transfer can be placed on the implant abutment prior to taking the impression in order to transfer the oral situation to the master cast.

However, insert.transfer is intended for single use only and cannot be re-sterilized. Hence a corresponding impression.transfer or set.cap must be used for subsequent impressions.

Packaging and sterility

ceramic implant dental implants are delivered in a sterile environment. The intact double-barrier sterile packaging protects the EO sterilised implant, including the accessories, from external influences, and if stored in a cool, dry place, guarantees sterility until the expiration date has elapsed. ceramic implant is intended for single use. Re-sterilisation is not possible. The contaminated implant should be disposed of together with residual waste.

After removing the ceramic implant from its sterile packaging, the instructions regarding asepsis should be followed. Open the packaging directly before implant insertion. Implants with damaged sterile packaging are not to be used. It is recommended to use the appropriate substitute product.

Further information

Detailed instructions for use and product information are available in the surgical guide and the prosthetic guide.

Explanation of symbols used

	Order number		Expiration date
	Batch designation		Store away from sunlight
	Sterile packaging, EO sterilised		0124 The CE mark confirms that the requirements of the Council Directive 93 / 42 EEC concerning Medical Devices are fulfilled
	Do not use in the case of damaged packaging		Attention – please refer to the enclosed documents
	Not for re-use		Manufacturer
	Date of manufacture		

Instructions for use for ceramic implant CI

vitaclinical is a business unit of VITA Zahnfabrik (referred to in the following as vitaclinical)
Manufacturer: VITA Zahnfabrik, Spitalgasse 3, 79713 Bad Säckingen, Germany

General

These instructions alone are not self-sufficient for the proper and safe use of the ceramic implant system. The product is only to be used by dental surgeons with a detailed knowledge of dental implantology and experience in the use of implant systems. It is absolutely essential to follow the instruction manual (surgical guide) and / or the instructions of an experienced operator. Training courses are offered by vitaclinical on demand.

Safety, liability and guarantee

The enclosed product is part of the vitaclinical implant system, and may only be used in conjunction with the corresponding original components and instruments as stipulated in the vitaclinical instructions and recommendations. The use of non-original, third-party components is detrimental to the performance of the implant, and leads to the loss of guarantee and replacement rights. Potential statutory liability claims for defects or product liability claims remain unaffected by this.

Instructions and advice regarding product application are transmitted in oral and written form, via electronic media or personal demonstration. This information, however, does not absolve the user from his or her personal responsibility to ensure that the product is suitable for its intended use, indications and procedures. The use of the product shall take place outside of the sphere of influence of vitaclinical, and hence falls into the user's area of responsibility.

vitaclinical does not accept any liability for direct or indirect consequential damages resulting from incorrect usage of the product. The liability of vitaclinical for potential defects is subject to the terms and conditions of sale and delivery of vitaclinical.

Copyright and trademark

vitaclinical or any other products and logos of vitaclinical mentioned here are trademarks or registered trademarks of VITA Zahnfabrik H. Rauter GmbH & Co. KG.

Product description

ceramic implant CI is a non-self-tapping ceramic dental implant. The screw thread is designed to generate a slight compression on the bone during insertion in order to create a high primary stability. The patented surface was especially developed for ceramic implant, cer.face 14, and helps accelerate osseointegration.

Indications

ceramic implant CI can be used for oral endosteal implants in the upper and lower jaw for single crowns and bridge abutments (for bridges with up to 3 units). It is suggested that prosthetic rehabilitation be carried out with all-ceramic single crowns or bridges (max.1 pontic). In the surgical guide you will find instructions on the recommended bone volume and quality and the distances between two implants, between implant and adjacent teeth.

Contraindications for ceramic implant

ceramic implant CI is not suitable for the following applications:

- Immediate implantation
- An edentulous jaw
- As anchorage for removable dentures
- In the case of insufficient bone volume
- In the case of insufficient bone quality resulting in insufficient primary implant stability

General relative contraindications

- Irradiated bone
- Diabetes mellitus
- Bruxism
- Parafunctions
- Anticoagulation drugs / haemorrhagic diathesis
- Crohn's disease
- Unsuitable anatomical bone situation
- Smoking
- Periodontitis
- Pathological alterations of the jaw and oral mucosa
- Pregnancy
- Insufficient oral hygiene